

# 10-Plex HCR™ Gold Spectral Imaging User Guide

This User Guide enables 10-Plex HCR™ Gold spectral imaging for robust imaging of 10 RNA and/or protein targets with 1-step quantitative HCR™ Gold signal amplification for all 10 targets simultaneously. 10-plex HCR™ Gold spectral imaging has the following key properties:

- No repeated imaging
  - Suitable even for whole-mounts and delicate specimens that are not amenable to repeated imaging
  - No repeated staining, imaging, stripping, or registration
- Standardized ingredients for robust performance
  - HCR™ Gold Amplifiers: orthogonal 10-set
  - Labels: optimized 10-set
  - Excitation wavelengths: optimized 10-set
  - Detection wavelengths: optimized 10-set
- Quantitative in all 10 channels
  - All 10 channels enable analog RNA or protein relative quantitation with subcellular resolution
  - All 10 channels enable digital RNA absolute quantitation with single-molecule resolution

## Contents

<b>10-Plex HCR™ Gold Spectral Imaging Workflow</b>	<b>2</b>
HCR™ Gold RNA-FISH/IF.....	2
AF Scan .....	2
Spectral Imaging.....	2
Linear Unmixing.....	3
<b>Frequently Asked Questions (FAQ)</b>	<b>4</b>
<b>HCR™ Technology Citation Notes</b>	<b>6</b>

## 10-Plex HCR™ Gold Spectral Imaging Workflow

Below is a general overview of the steps involved in 10-Plex HCR™ Gold spectral imaging. See Supplementary Section S2.1 of [Schulte et al., 2024](#) for a detailed protocol.

### HCR™ Gold RNA-FISH/IF

Follow the protocols in the appropriate HCR™ Gold User Guides:

- [HCR™ Gold RNA-FISH User Guide](#)
- [HCR™ Gold IF User Guide](#)
- [HCR™ Gold RNA-FISH/IF User Guide](#)

to perform HCR™ Gold RNA-FISH and/or IF including:

- Sample preparation
- Antibody binding (for protein targets)
- Probe hybridization (for RNA targets)
- Amplification

on 12 sample types:

- 10-plex sample or multiple replicate samples (1 label for each of 10 targets): use all 10 HCR™ Gold Probes and Amplifiers.
- 1-plex reference sample for each of 10 targets: use the corresponding HCR™ Gold Probe and Amplifier for a given target.
- Two unlabeled autofluorescence (AF) samples: omit all HCR™ Gold Probes and Amplifiers.

**NOTE:** When assigning targets to Ch1-Ch10, we recommend assigning higher-expression targets to channels that employ labels with emissions spectra strongly overlapping sample autofluorescence (e.g., Ch1-Ch3 for whole-mount zebrafish embryos and fresh-frozen mouse brain sections) or with emissions spectra for which the microscope detectors are less sensitive (e.g., Ch10 in the near-IR wavelength range).

### AF Scan

Use one unlabeled sample to perform an excitation-emission scan to determine the maximal AF excitation wavelength ( $\lambda_{AF}$ ), which in turn determines the optimized detection wavelengths for AF (see Table 1).

**NOTE:** approximately 30–45 sec of microscope time

### Spectral Imaging

Spectrally image 12 sample types using 11 excitation wavelengths in succession (one optimized for each label and one optimized for AF; see Table 1):

- 10-plex sample.
- 1-plex samples (obtain reference sample for each label).
- Unlabeled sample (obtain reference spectrum for AF).

**NOTE:** approximately 30–45 sec of microscope time per optical section

## Linear Unmixing

Use the 11 reference spectra (one per label and one for AF) to linearly unmix the 10-plex image to obtain 11 unmixed channels (one per label and one for AF).

**NOTE:** *approximately 5-10 sec of computational time*

**Table 1.** Excitation and detection wavelengths for 10-plex HCR™ Gold spectral imaging.

Channel	Excitation wavelength (nm)	Detector wavelengths (nm)	Label
1	405	410–450	405
2	440	450–520	425
3	488	493–533	488
4	518	523–583	514
5	557	566–620	546
6	590	600–660	594
7	629	640–680	633
8	686	696–723	700
9	755	765–795	750
10	790	815–850	800
AF	$\lambda_{AF}$	$(\lambda_{AF} + 6) - (\lambda_{AF} + 171)$	–

## Frequently Asked Questions (FAQ)

### What comes in an HCR™ Gold RNA-FISH kit?

- HCR™ HiFi Probe
- HCR™ HiFi Probe Hybridization Buffer
- HCR™ HiFi Probe Wash Buffer
- HCR™ Gold Amplifier
- HCR™ Gold Amplifier Buffer
- HCR™ Gold Amplifier Wash Buffer

### What comes in an HCR™ Gold IF kit?

- HCR™ 2° Antibody Probe
- HCR™ Antibody Buffer
- HCR™ Gold Amplifier
- HCR™ Gold Amplifier Buffer
- HCR™ Gold Amplifier Wash Buffer

### Can I order a subset of these components?

- Yes. You can order kits containing any subset of the above components.

### What HCR™ Gold Amplifiers are available for multiplexing?

- For multiplex experiments, choose a different HCR™ Gold Amplifier (X1, X2, ..., X10) and label (405, 425, 488, 514, 546, 594, 633, 647, 700, 750, 800) for each target that will be imaged in the same sample (for example, amplifier X3-647 for target 1, amplifier X2-546 for target 2, ...).
- We recommend labels (488, 546, 647, 750) for robust 4-plex bandpass imaging.
- We recommend labels (405, 425, 488, 514, 546, 594, 633, 700, 750, 800) for robust 10-plex spectral imaging.
- For flexibility, amplifiers X1, X2, and X3 are available with multiple label options (488, 546, 647), while for expanded multiplexing, amplifiers X4–X10 are each available with one label option.
- We recommend using a longer-wavelength label (e.g., 647 or 750) for targets that are more difficult to detect (e.g., low-expression target and/or short target sequence), as autofluorescence tends to be higher in shorter-wavelength channels (e.g., 405 or 488).

### What do I order for a multiplex HCR™ Gold RNA-FISH/IF experiment?

- Order one HCR™ Gold RNA-FISH kit for each target RNA
- Order one HCR™ Gold IF kit for each target protein

- Example 4-Plex Experiment (2 RNAs, 2 proteins):
  - HCR™ Gold RNA-FISH kit for target RNA1:
    - HCR™ HiFi Probe: target RNA1 for use with amplifier X3
    - HCR™ HiFi Probe Hybridization Buffer (for use with all RNA-FISH kits)
    - HCR™ HiFi Probe Wash Buffer (for use with all RNA-FISH kits)
    - HCR™ Gold Amplifier: X3 with label 647
    - HCR™ Gold Amplifier Buffer (for use with all RNA-FISH and IF kits)
    - HCR™ Gold Amplifier Wash Buffer (for use with all RNA-FISH and IF kits)
  - HCR™ Gold RNA-FISH kit for target RNA2:
    - HCR™ HiFi Probe: target RNA2 for use with amplifier X4
    - HCR™ Gold Amplifier: X4 with label 750
  - HCR™ Gold IF kit for target Protein1:
    - 1° antibody: Rabbit Anti-Protein1 (your own 1° antibody)
    - HCR™ 2° Antibody Probe: Donkey Anti-Rabbit for use with amplifier X1
    - HCR™ Antibody Buffer (for use with all IF kits)
    - HCR™ Gold Amplifier: X1 with label 488
  - HCR™ Gold IF kit for target Protein2:
    - 1° antibody: Mouse Anti-Protein2 (your own 1° antibody)
    - HCR™ 2° Antibody Probe: Donkey Anti-Mouse for use with amplifier X2
    - HCR™ Gold Amplifier: X2 with label 546

**Are HCR™ Gold Amplifiers interchangeable between HCR™ Gold RNA-FISH kits and HCR™ Gold IF kits?**

- Yes! The HCR™ Gold Amplifiers (488, 546, 647) are interchangeable between HCR™ Gold RNA-FISH and HCR™ Gold IF kits. Likewise, HCR™ Gold Amplifier Buffer and HCR™ Gold Amplifier Wash Buffer are interchangeable between HCR™ Gold RNA-FISH and HCR™ Gold IF kits.

## HCR™ Technology Citation Notes

For citation, please select from the list below as appropriate for your application:

- **10-Plex HCR™ Spectral Imaging**  
HCR™ RNA-FISH/IF enables quantitative high-resolution imaging of 10 RNA and/or protein targets with 1-step HCR™ signal amplification for all targets simultaneously. The method is suitable even for whole-mounts and delicate samples as it requires no repeated staining, imaging, registration, or stripping ([Schulte et al., 2024](#)).
- **HCR™ RNA-FISH/IF**  
HCR™ RNA-FISH/IF enables a unified approach to multiplex, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) and protein immunofluorescence (IF), with quantitative 1-step enzyme-free signal amplification performed for all RNA and protein targets simultaneously ([Schwarzkopf et al., 2021](#)).
- **HCR™ IF**  
HCR™ IF enables multiplex, quantitative, high-resolution protein immunofluorescence (IF) in highly autofluorescent samples (e.g., FFPE brain tissue sections) ([Schwarzkopf et al., 2021](#)).
- **HCR™ RNA-FISH**
  - Third-generation HCR™ RNA-FISH (v3.0) enables multiplex, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) with automatic background suppression throughout the protocol for dramatically enhanced performance (signal-to-background, subcellular quantitative RNA imaging precision, single-molecule quantitative RNA imaging fidelity) and ease-of-use (no probe set optimization for new targets and organisms) ([Choi et al., 2018](#)).
  - Second-generation HCR™ RNA-FISH (v2.0) using DNA HCR™ Probes and DNA HCR™ Amplifiers: 10× increase in signal, 10× reduction in cost, dramatic increase in reagent durability ([Choi et al., 2014](#)).
  - First-generation HCR™ RNA-FISH (v1.0) using RNA HCR™ Probes and RNA HCR™ Amplifiers: multiplex mRNA imaging in whole-mount vertebrate embryos with simultaneous signal amplification for up to 5 target mRNAs ([Choi et al., 2010](#)).
- **Subcellular Quantitative RNA and Protein Imaging**  
HCR™ RNA-FISH enables analog relative quantitation of RNA and/or protein targets with subcellular resolution in the anatomical context of thick autofluorescent samples (e.g., whole-mount vertebrate embryos) ([Trivedi et al., 2018](#), [Choi et al., 2018](#), [Schwarzkopf et al., 2021](#)).
- **Single-Molecule Quantitative RNA Imaging**  
HCR™ RNA-FISH enables digital RNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) ([Shah et al., 2016](#), [Choi et al., 2018](#)).
- **Read-Out/Read-In Analysis Framework**  
The read-out/read-in analysis framework enables bidirectional quantitative discovery in an anatomical context ([Trivedi et al., 2018](#)).
- **Protocols in Diverse Sample Types**  
Protocols for HCR™ RNA-FISH and/or IF in diverse sample types are adapted from the zoo paper ([Choi et al., 2016](#)):

- bacteria in suspension
  - FFPE human tissue sections
  - generic sample in solution
  - generic sample on a slide
  - mammalian cells on a slide
  - mammalian cells in suspension
  - whole-mount chicken embryos
  - whole-mount fruit fly embryos
  - whole-mount mouse embryos
  - whole-mount nematode larvae
  - whole-mount sea urchin embryos
  - whole-mount zebrafish embryos and larvae
- **HCR™ RNA Flow Cytometry**  
HCR™ RNA Flow Cytometry enables analog RNA relative quantitation for high-throughput expression profiling of mammalian cells and bacteria without the need to engineer reporter lines ([Choi et al., 2018](#)).
  - **HCR™ Northern Blots**  
HCR™ Northern Blots enable simultaneous quantification of RNA target size and abundance with automatic background suppression throughout the protocol ([Schwarzkopf & Pierce, 2016](#)).
  - **HCR™ Amplifiers**  
HCR™ Amplifiers enable multiplex, quantitative, 1-step, isothermal, enzyme-free signal amplification in diverse technological settings ([Dirks & Pierce, 2004](#)).