

## HCR v3.0 protocol for whole-mount fruit fly embryos (*Drosophila melanogaster*)

This protocol has not been optimized for all stages and should only be used as a template.



Technical support: [support@molecularinstruments.com](mailto:support@molecularinstruments.com)

Safety data sheets (SDS): [www.molecularinstruments.com/safety](http://www.molecularinstruments.com/safety)

Molecular Instruments products are protected by and for use under patents: [www.molecularinstruments.com/patents](http://www.molecularinstruments.com/patents)

### Storage conditions

Store HCR probe sets, HCR amplifiers, HCR probe hybridization buffer, and HCR probe wash buffer at -20 °C. Store HCR amplification buffer at 4 °C. On the bench top, keep stock solutions on ice. Make sure all solutions are well mixed before use.



## Preparation of fixed whole-mount fruit fly embryos

1. Collect fly embryos and incubate with yeast paste (food source) until they reach stage 4–6 (approximately 3 h).
2. Rinse embryos into a mesh basket and wash excess yeast paste using DI H<sub>2</sub>O.
3. Immerse the mesh basket into 100% bleach and wash for 2 min to dechorionate embryos.
4. Rinse the basket with DI H<sub>2</sub>O.
5. Transfer embryos to a scintillation vial containing 8 mL of 4.5% formaldehyde fixation solution.  
*CAUTION: use formaldehyde with extreme care as it is a hazardous material.*
6. Gently rock embryos in scintillation vial for 25 min on an orbital rocker.  
*NOTE: Rocking does not have to be very vigorous. However, embryos must move continuously within the heptane-fixative interphase to be exposed uniformly to the solvent and fixative.*
7. Remove the bottom liquid phase in the vial.
8. Add 8 mL of methanol (MeOH) and shake the vial hard for 1 min. Devitellinized embryos will sink to the bottom of the vial.
9. Remove all liquid and rinse 2 times in 1 mL of MeOH to remove debris.  
*NOTE: Embryos can be stored in 1 mL of MeOH at -20 °C before use.*
10. For each sample, transfer 50 μL of embryos (using a cut pipet tip) to a 1.5 mL tube.
11. Rinse embryos 4 times by quickly adding 1 mL of ethanol (EtOH), inverting the tube, and aspirating the supernatant.
12. Add 500 μL of EtOH and 250 μL of xylene and invert the tube.  
*CAUTION: use xylene with care as it is a hazardous material.*
13. Add an additional 250 μL of xylene and invert the tube.
14. Add another 250 μL of xylene again and invert the tube.  
*NOTE: The tube should now contain 500 μL of EtOH and 750 μL of xylene.*
15. Rock the tube at room temperature for 45–60 min. Longer rocking time is acceptable.
16. Aspirate the supernatant.
17. Rinse embryos once and wash 3 × 5 min with EtOH.  
*NOTE: All washes have a volume of 1 mL unless specified. All washes before pre-hybridization (step 1) are performed with rocking.*
18. Rinse embryos once and wash 2 × 5 min with MeOH.

19. Wash with 50% MeOH / 50% PBST for 5 min.
20. Wash 1 × 10 min and 2 × 5 min with PBST.
21. Rock embryos in 1 mL of 4 μg/mL proteinase K solution at room temperature for 7 min.  
*NOTE: Proteinase K concentration and treatment time should be reoptimized for each batch of proteinase K, or for samples at a different developmental stage.*
22. Rinse embryos 2 times and wash 2 × 5 min with PBST.
23. Rock embryos in 4% formaldehyde post-fixation solution at room temperature for 25 min.
24. Rinse embryos and wash 5 × 5 min with PBST.

## Buffer recipes for sample preparation

### 4.5% formaldehyde fixation solution

4.5% formaldehyde  
0.5× PBS  
25 mM EGTA  
50% heptane

### For 8 mL of solution

975 μL of 37% formaldehyde  
400 μL of 10× PBS  
76 mg of EGTA  
4 mL of heptane  
Fill up to 8 mL with ultrapure H<sub>2</sub>O

### Proteinase K solution

4 μg/mL proteinase K

### For 2 mL of solution

0.4 μL of 20 mg/mL proteinase K  
Fill up to 2 mL with PBST

### 4% formaldehyde post-fixation solution

4% formaldehyde

### For 2 mL of solution

216 μL of 37% formaldehyde  
Fill up to 2 mL with PBST

### PBST

1× PBS  
0.1% Tween 20

### For 50 mL of solution

5 mL of 10× PBS  
500 μL of 10% Tween 20  
Fill up to 50 mL with ultrapure H<sub>2</sub>O

*NOTE: avoid using calcium chloride and magnesium chloride in PBS as this leads to increased autofluorescence in the embryos.*

## Multiplexed HCR v3.0 protocol

### Detection stage

1. Pre-hybridize embryos in 200  $\mu\text{L}$  of probe hybridization buffer for 30 min at 37  $^{\circ}\text{C}$ .  
*CAUTION: probe hybridization buffer contains formamide, a hazardous material.*
2. Prepare probe solution by adding 0.8 pmol of each probe set (e.g. 0.8  $\mu\text{L}$  of 1  $\mu\text{M}$  stock) to 200  $\mu\text{L}$  of probe hybridization buffer at 37  $^{\circ}\text{C}$ .  
*NOTE: For dHCR imaging, use higher concentration of probe to improve probe hybridization efficiency.*
3. Remove the pre-hybridization solution and add the probe solution.
4. Incubate embryos overnight (12–16 h) at 37  $^{\circ}\text{C}$ .
5. Remove excess probes by washing embryos 4  $\times$  15 min with 1 mL of probe wash buffer at 37  $^{\circ}\text{C}$ .  
*CAUTION: probe wash buffer contains formamide, a hazardous material.*  
*NOTE: pre-heat probe wash buffer to 37  $^{\circ}\text{C}$  before use.*
6. Wash samples 2  $\times$  5 min with 1 mL of 5 $\times$  SSCT at room temperature.

### Amplification stage

1. Pre-amplify embryos with 1 mL of amplification buffer for 10 min at room temperature.  
*NOTE: equilibrate amplification buffer to room temperature before use.*
2. Separately prepare 6 pmol of hairpin h1 and 6 pmol of hairpin h2 by snap cooling 2  $\mu\text{L}$  of 3  $\mu\text{M}$  stock (heat at 95  $^{\circ}\text{C}$  for 90 seconds and cool to room temperature in a dark drawer for 30 min).  
*NOTE: HCR hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. h1 and h2 should be snap cooled in separate tubes.*
3. Prepare hairpin solution by adding snap-cooled h1 hairpins and snap-cooled h2 hairpins to 100  $\mu\text{L}$  of amplification buffer at room temperature.
4. Remove the pre-amplification solution and add the hairpin solution.
5. Incubate the embryos overnight (12–16 h) in the dark at room temperature.  
*NOTE: For dHCR imaging, amplify for a shorter period of time to ensure single-molecule dots are diffraction-limited.*
6. Remove excess hairpins by washing with 1 mL of 5 $\times$  SSCT at room temperature:
  - (a) 2  $\times$  5 min
  - (b) 2  $\times$  30 min
  - (c) 1  $\times$  5 min
7. Samples can be stored at 4  $^{\circ}\text{C}$  protected from light before microscopy.

## Citation Notes

For citation, please select from the list below as appropriate for your application:

- **HCR v3.0**  
Automatic background suppression for dramatically enhanced performance (signal-to-background, qHCR precision, dHCR fidelity) and ease-of-use (no probe set optimization for new targets and organisms) (Choi *et al.*, 2018). Quantitative analysis modes:
  - **qHCR imaging**: analog mRNA relative quantitation with subcellular resolution in the anatomical context of thick autofluorescent samples.
  - **qHCR flow cytometry**: analog mRNA relative quantitation for high-throughput expression profiling of mammalian and bacterial cells.
  - **dHCR imaging**: digital mRNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples.

Protocols for v3.0 in diverse organisms are adapted from the Zoo paper.
- **qHCR imaging**  
mRNA relative quantitation with subcellular resolution in the anatomical context of whole-mount vertebrate embryos; read-out/read-in analysis framework (Trivedi *et al.*, 2018).
- **Zoo paper**  
Protocols for multiplexed mRNA imaging in diverse sample types (Choi *et al.*, 2016):
 

bacteria in suspension	FFPE human tissue sections
generic sample in suspension	generic sample on slide
whole-mount chicken embryos	whole-mount fruit fly embryos
whole-mount mouse embryos	whole-mount sea urchin embryos
whole-mount worm larvae	whole-mount zebrafish embryos and larvae
- **dHCR imaging**  
Single-molecule mRNA imaging in thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) (Shah *et al.*, 2016).
- **qHCR northern blot**  
Simultaneous quantification of RNA target size and abundance for up to 5 target RNAs (Schwarzkopf & Pierce, 2016).
- **HCR v2.0**  
2nd generation in situ HCR technology (v2.0) using DNA HCR probes and DNA HCR amplifiers: 10× increase in signal, 10× reduction in cost, dramatic increase in reagent durability (Choi *et al.*, 2014).
- **HCR v1.0**  
1st generation in situ HCR technology (v1.0) using RNA HCR probes and RNA HCR amplifiers: multiplexed mRNA imaging in whole-mount vertebrate embryos with simultaneous signal amplification for up to 5 target mRNAs (Choi *et al.*, 2010).
- **HCR mechanism**  
Isothermal enzyme-free molecular signal amplification (Dirks & Pierce, 2004).

## References

- Choi, H. M. T., Chang, J. Y., Trinh, L. A., Padilla, J. E., Fraser, S. E., & Pierce, N. A. (2010). Programmable in situ amplification for multiplexed imaging of mRNA expression. *Nat Biotechnol*, **28**(11), 1208–12.
- Choi, H. M. T., Beck, V. A., & Pierce, N. A. (2014). Next-generation in situ hybridization chain reaction: higher gain, lower cost, greater durability. *ACS Nano*, **8**(5), 4284–4294.
- Choi, H. M. T., Calvert, C. R., Husain, N., Huss, D., Barsi, J. C., Deverman, B. E., Hunter, R. C., Kato, M., Lee, S. M., Abelin, A. C. T., Rosenthal, A. Z., Akbari, O. S., Li, Y., Hay, B. A., Sternberg, P. W., Patterson, P. H., Davidson, E. H., Mazmanian, S. K., Prober, D. A., van de Rijn, M., Leadbetter, J. R., Newman, D. K., Readhead, C., Bronner, M. E., Wold, B., Lansford, R., Sauka-Spengler, T., Fraser, S. E., & Pierce, N. A. (2016). Mapping a multiplexed zoo of mRNA expression. *Development*, **143**, 3632–3637.
- Choi, H. M. T., Schwarzkopf, M., Fornace, M. E., Acharya, A., Artavanis, G., Stegmaier, J., Cunha, A., & Pierce, N. A. (2018). Third-generation in situ hybridization chain reaction: multiplexed, quantitative, sensitive, versatile, robust. *Development*, **145**, dev165753.
- Dirks, R. M., & Pierce, N. A. (2004). Triggered amplification by hybridization chain reaction. *Proc Natl Acad Sci USA*, **101**(43), 15275–15278.
- Schwarzkopf, M., & Pierce, N. A. (2016). Multiplexed miRNA northern blots via hybridization chain reaction. *Nucleic Acids Res*, **44**(15), e129.
- Shah, S., Lubeck, E., Schwarzkopf, M., He, T.-F., Greenbaum, A., Sohn, C. H., Lignell, A., Choi, H. M. T., Gradi-naru, V., Pierce, N. A., & Cai, L. (2016). Single-molecule RNA detection at depth via hybridization chain reaction and tissue hydrogel embedding and clearing. *Development*, **143**, 2862–2867.
- Trivedi, V., Choi, H. M. T., Fraser, S. E., & Pierce, N. A. (2018). Multidimensional quantitative analysis of mRNA expression within intact vertebrate embryos. *Development*, **145**, dev.156869.