

## HCR v3.0 protocol for fresh frozen or fixed frozen tissue sections

This protocol has not been validated for all sample types yet and should only be used as a template.



Technical support: [support@molecularinstruments.com](mailto:support@molecularinstruments.com)

Safety data sheets (SDS): [www.molecularinstruments.com/safety](http://www.molecularinstruments.com/safety)

Molecular Instruments products are protected by and for use under patents: [www.molecularinstruments.com/patents](http://www.molecularinstruments.com/patents)

### Storage conditions

Store HCR probe sets, HCR amplifiers, HCR probe hybridization buffer, and HCR probe wash buffer at -20 °C. Store HCR amplification buffer at 4 °C. On the bench top, keep stock solutions on ice. Make sure all solutions are well mixed before use.



## Sample preparation protocol

1. Remove frozen sections on slide from -80 °C.
2. Fix tissues by immersing slides in ice-cold 4% paraformaldehyde (PFA) for 15 min at 4 °C.  
*CAUTION: use PFA with extreme care as it is a hazardous material.*  
*NOTE: use fresh PFA and cool to 4 °C before use to avoid increased autofluorescence.*  
*NOTE: Each 50 mL tube can fit two outward-facing slides. A volume of 30 mL is sufficient to immerse sections in the tube. If desired, a larger number of slides can be processed together using a Coplin jar.*
3. Immerse slides in 50% EtOH for 5 min at room temperature.
4. Immerse slides in 70% EtOH for 5 min at room temperature.
5. Immerse slides in 100% EtOH for 5 min at room temperature.
6. Immerse slides in another fresh 100% EtOH for 5 min at room temperature.
7. Remove slides and allow them to air dry for 5 min at room temperature.
8. Introduce 200  $\mu$ L of 10  $\mu$ g/mL of proteinase K solution on top of sample and incubate slides in a humidified chamber for 10 min at 37 °C.  
*NOTE: Proteinase K concentration and treatment time should be reoptimized for each batch of proteinase K.*
9. Wash slides by immersing in 1 $\times$  PBS.  
*NOTE: avoid using calcium chloride and magnesium chloride in PBS as this leads to increased autofluorescence in the tissue.*
10. Repeat with fresh 1 $\times$  PBS.
11. Proceed to HCR assay.

## Buffer recipes for sample preparation

### 4% paraformaldehyde (PFA)

4% PFA

1× PBS

### For 30 mL of solution

7.5 mL of 16% PFA solution

3 mL of 10× PBS

Fill up to 30 mL with water

### Proteinase K solution

10  $\mu$ g/mL proteinase K

### For 1 mL of solution

0.5  $\mu$ L of 20 mg/mL proteinase K

Fill up to 1 mL with 1× PBS

*NOTE: avoid using calcium chloride and magnesium chloride in PBS as this leads to increased autofluorescence in the samples.*

## Multiplexed HCR v3.0 protocol

### Detection stage

1. Pre-warm a humidified chamber to 37 °C.
2. Dry slide by blotting edges on a Kimwipe.
3. Add 200  $\mu$ L of probe hybridization buffer on top of the sample.  
*CAUTION: probe hybridization buffer contains formamide, a hazardous material.*
4. Pre-hybridize for 10 min inside the humidified chamber.
5. Prepare probe solution by adding 0.4 pmol of each probe set (e.g. 0.4  $\mu$ L of 1  $\mu$ M stock) to 100  $\mu$ L of probe hybridization buffer at 37 °C.  
*NOTE: For dHCR imaging, use higher concentration of probe to improve probe hybridization efficiency.*
6. Remove the pre-hybridization solution and drain excess buffer on slide by blotting edges on a Kimwipe.
7. Add 50–100  $\mu$ L of the probe solution on top of the sample.  
*NOTE: Amount of probe solution depends on the size of the coverslip.*
8. Place a coverslip on the sample and incubate overnight (12–16 h) in the 37 °C humidified chamber.
9. Immerse slide in probe wash buffer at 37 °C to float off coverslip.  
*CAUTION: probe wash buffer contains formamide, a hazardous material.*
10. Remove excess probes by incubating slide at 37 °C in:
  - (a) 75% of probe wash buffer / 25% 5 $\times$  SSCT for 15 min
  - (b) 50% of probe wash buffer / 50% 5 $\times$  SSCT for 15 min
  - (c) 25% of probe wash buffer / 75% 5 $\times$  SSCT for 15 min
  - (d) 100% 5 $\times$  SSCT for 15 min*NOTE: Wash solutions should be pre-heated to 37 °C before use.*
11. Immerse slide in 5 $\times$  SSCT for 5 min at room temperature.

## Amplification stage

1. Dry slide by blotting edges on a Kimwipe.
2. Add 200  $\mu\text{L}$  of amplification buffer on top of the sample and pre-amplify in a humidified chamber for 30 min at room temperature.
3. Separately prepare 6 pmol of hairpin h1 and 6 pmol of hairpin h2 by snap cooling 2  $\mu\text{L}$  of 3  $\mu\text{M}$  stock (heat at 95  $^{\circ}\text{C}$  for 90 seconds and cool to room temperature in a dark drawer for 30 min).  
*NOTE: HCR hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. h1 and h2 should be snap cooled in separate tubes.*
4. Prepare hairpin solution by adding snap-cooled h1 hairpins and snap-cooled h2 hairpins to 100  $\mu\text{L}$  of amplification buffer at room temperature.
5. Remove the pre-amplification solution and drain excess buffer on slide by blotting edges on a Kimwipe.
6. Add 50–100  $\mu\text{L}$  of the hairpin solution on top of the sample.  
*NOTE: Amount of hairpin solution depends on the size of the coverslip.*
7. Place a coverslip on the sample and incubate overnight (12–16 h) in a dark humidified chamber at room temperature.  
*NOTE: For dHCR imaging, amplify for a shorter period of time to ensure single-molecule dots are diffraction-limited.*
8. Immerse slide in  $5 \times$  SSCT at room temperature to float off coverslip.
9. Remove excess hairpins by incubating slide in  $5 \times$  SSCT at room temperature for:
  - (a)  $2 \times 30$  min
  - (b)  $1 \times 5$  min
10. Dry slide by blotting edges on a Kimwipe.
11. Add 50–100  $\mu\text{L}$  of antifade mounting reagent on top of the sample.
12. Place a coverslip on top for microscopy.

## Citation Notes

For citation, please select from the list below as appropriate for your application:

- **HCR v3.0**

Automatic background suppression for dramatically enhanced performance (signal-to-background, qHCR precision, dHCR fidelity) and ease-of-use (no probe set optimization for new targets and organisms) (Choi *et al.*, 2018). Quantitative analysis modes:

**qHCR imaging:** analog mRNA relative quantitation with subcellular resolution in the anatomical context of thick autofluorescent samples.

**qHCR flow cytometry:** analog mRNA relative quantitation for high-throughput expression profiling of mammalian and bacterial cells.

**dHCR imaging:** digital mRNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples.

Protocols for v3.0 in diverse organisms are adapted from the Zoo paper.

- **qHCR imaging**

mRNA relative quantitation with subcellular resolution in the anatomical context of whole-mount vertebrate embryos; read-out/read-in analysis framework (Trivedi *et al.*, 2018).

- **Zoo paper**

Protocols for multiplexed mRNA imaging in diverse sample types (Choi *et al.*, 2016):

bacteria in suspension	FFPE human tissue sections
generic sample in suspension	generic sample on slide
whole-mount chicken embryos	whole-mount fruit fly embryos
whole-mount mouse embryos	whole-mount sea urchin embryos
whole-mount worm larvae	whole-mount zebrafish embryos and larvae

- **dHCR imaging**

Single-molecule mRNA imaging in thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) (Shah *et al.*, 2016).

- **qHCR northern blot**

Simultaneous quantification of RNA target size and abundance for up to 5 target RNAs (Schwarzkopf & Pierce, 2016).

- **HCR v2.0**

2nd generation in situ HCR technology (v2.0) using DNA HCR probes and DNA HCR amplifiers: 10× increase in signal, 10× reduction in cost, dramatic increase in reagent durability (Choi *et al.*, 2014).

- **HCR v1.0**

1st generation in situ HCR technology (v1.0) using RNA HCR probes and RNA HCR amplifiers: multiplexed mRNA imaging in whole-mount vertebrate embryos with simultaneous signal amplification for up to 5 target mRNAs (Choi *et al.*, 2010).

- **HCR mechanism**

Isothermal enzyme-free molecular signal amplification (Dirks & Pierce, 2004).

## References

- Choi, H. M. T., Chang, J. Y., Trinh, L. A., Padilla, J. E., Fraser, S. E., & Pierce, N. A. (2010). Programmable in situ amplification for multiplexed imaging of mRNA expression. *Nat Biotechnol*, **28**(11), 1208–12.
- Choi, H. M. T., Beck, V. A., & Pierce, N. A. (2014). Next-generation in situ hybridization chain reaction: higher gain, lower cost, greater durability. *ACS Nano*, **8**(5), 4284–4294.
- Choi, H. M. T., Calvert, C. R., Husain, N., Huss, D., Barsi, J. C., Deverman, B. E., Hunter, R. C., Kato, M., Lee, S. M., Abelin, A. C. T., Rosenthal, A. Z., Akbari, O. S., Li, Y., Hay, B. A., Sternberg, P. W., Patterson, P. H., Davidson, E. H., Mazmanian, S. K., Prober, D. A., van de Rijn, M., Leadbetter, J. R., Newman, D. K., Readhead, C., Bronner, M. E., Wold, B., Lansford, R., Sauka-Spengler, T., Fraser, S. E., & Pierce, N. A. (2016). Mapping a multiplexed zoo of mRNA expression. *Development*, **143**, 3632–3637.
- Choi, H. M. T., Schwarzkopf, M., Fornace, M. E., Acharya, A., Artavanis, G., Stegmaier, J., Cunha, A., & Pierce, N. A. (2018). Third-generation in situ hybridization chain reaction: multiplexed, quantitative, sensitive, versatile, robust. *Development*, **145**, dev165753.
- Dirks, R. M., & Pierce, N. A. (2004). Triggered amplification by hybridization chain reaction. *Proc Natl Acad Sci USA*, **101**(43), 15275–15278.
- Schwarzkopf, M., & Pierce, N. A. (2016). Multiplexed miRNA northern blots via hybridization chain reaction. *Nucleic Acids Res*, **44**(15), e129.
- Shah, S., Lubeck, E., Schwarzkopf, M., He, T.-F., Greenbaum, A., Sohn, C. H., Lignell, A., Choi, H. M. T., Gradi-naru, V., Pierce, N. A., & Cai, L. (2016). Single-molecule RNA detection at depth via hybridization chain reaction and tissue hydrogel embedding and clearing. *Development*, **143**, 2862–2867.
- Trivedi, V., Choi, H. M. T., Fraser, S. E., & Pierce, N. A. (2018). Multidimensional quantitative analysis of mRNA expression within intact vertebrate embryos. *Development*, **145**, dev.156869.