

HCR™ northern blot (v3.0) protocol

Technical Support

support@molecularinstruments.com

Safety Data Sheets (SDS)

www.molecularinstruments.com/safety-v3

Patents

www.molecularinstruments.com/patents

Ordering for Multiplex Experiment

Order one HCR™ RNA-FISH (v3.0) kit per target RNA

Example 2-Plex Experiment

- HCR™ RNA-FISH (v3.0) kit for target mRNA1
 - HCR™ Probe (v3.0): target mRNA1 for use with amplifier B1
 - HCR™ Amplifier (v3.0): B1-647
 - HCR™ RNA-FISH Buffers (v3.0): HCR™ Probe Hybridization Buffer (v3.0), HCR™ Probe Wash Buffer (v3.0), HCR™ Amplifier Buffer (v3.0) (for use with all kits)
- HCR™ RNA-FISH (v3.0) kit for target mRNA2
 - HCR™ Probe (v3.0): target mRNA2 for use with amplifier B2
 - HCR™ Amplifier (v3.0): B2-488

Storage conditions

- Store HCR™ Probes (v3.0), HCR™ Amplifiers (v3.0), HCR™ Probe Hybridization Buffer (v3.0), and HCR™ Probe Wash Buffer (v3.0) at -20 °C.
- Store HCR™ Amplifier Buffer (v3.0) at 4 °C.
- On the bench top, keep stock solutions on ice.
- Make sure all solutions are well mixed before use.

Denaturing agarose gel electrophoresis

1. Prepare a 1% denaturing agarose gel.
NOTE: The gel should be ≤ 6 mm for efficient RNA transfer.
NOTE: We recommend using NorthernMax™ denaturing gel buffer for preparing denaturing agarose gels.
2. Mix RNA samples in a 1:1 ratio with formamide and heat to 65 °C for 15 min prior to gel loading.
CAUTION: formamide is a hazardous material.
3. Run RNA samples (typically 1–20 μ g) at 5V per cm in 1 \times NorthernMax™ running buffer until your targets of interest are well-resolved.
NOTE: The distance measured is the distance between the electrodes in the electrophoresis chamber.
NOTE: For targets with low or unknown expression levels we recommend starting with 10–20 μ g of total RNA.

RNA transfer and crosslinking

1. Wash the gel 4 \times 5 min, in ultrapure H₂O.
2. Wet membrane in ultrapure H₂O and equilibrate for 5 min in 20 \times SSC.
NOTE: A positively charged nylon membrane is recommended.
3. Perform a capillary transfer of RNA onto the membrane in 20 \times SSC.
NOTE: We recommend using the TurboBlotter Rapid Downward Transfer System. Alternative transfer methods may be used.
4. Wash the membrane for 5 min in 2 \times SSC.
5. Crosslink the RNA to the membrane by incubating at 80 °C for 2 h.
NOTE: Alternative RNA crosslinking methods may be used.

Multiplexed HCR™ northern blot (v3.0) protocol

Detection stage

1. Let V denote the volume of buffer to be used for blots, calculated as:

- (a) 1 mL per 10 cm² of membrane
- (b) rounded up to the nearest 0.5 mL
- (c) no less than 1.5 mL

For example, $V = 3.5$ mL for a membrane that is 32 cm² and $V = 1.5$ mL for a membrane that is 9 cm².

2. Pre-heat HCR™ Probe Hybridization Buffer (v3.0) to 37 °C.
CAUTION: HCR™ Probe Hybridization Buffer (v3.0) contains formamide, a hazardous material.
3. Pre-hybridize blot in ($V - 0.5$ mL) of HCR™ Probe Hybridization Buffer (v3.0) in a hybridization bottle for 30–60 min at 37 °C in a rolling hybridization incubator.
4. Prepare probe solution by adding each HCR™ Probe (v3.0) into a total volume of 0.5 mL of HCR™ Probe Hybridization Buffer (v3.0) at 37 °C such that each probe will be at a final concentration of 5 nM in volume V .
5. Add the probe solution to the pre-hybridization solution and incubate blot overnight at 37 °C in a rolling hybridization incubator.
6. Remove excess probes by washing four times using pre-heated 37 °C wash solutions (volume $2V$ per wash) :
 - (a) two low-stringency washes ($2 \times$ SSC, 0.1% SDS; 5 min at 37 °C)
 - (b) two high-stringency washes ($0.2 \times$ SSC, 0.1% SDS; 15 min at 37 °C)

Amplification stage

1. Pre-heat HCR™ Amplifier Buffer (v3.0) to 37 °C.
2. Pre-amplify blot in ($V - 0.5$ mL) of HCR™ Amplifier Buffer (v3.0) in a hybridization bottle for 30–60 min at 37 °C in a rolling hybridization incubator.
3. Separately prepare hairpin h1 and hairpin h2 by snap cooling (heat at 95 °C for 90 seconds and cool to room temperature in a dark drawer for 30 min).
NOTE: Hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. Snap cool an amount of each hairpin corresponding to 30 nM final concentration in a final volume V . h1 and h2 should be snap cooled in separate tubes.
4. Prepare hairpin solution by adding all snap-cooled hairpins into a final volume of 0.5 mL HCR™ Amplifier Buffer (v3.0) at 37 °C.
5. Add the hairpin solution to the pre-amplification solution and incubate blot for 4 h at 37 °C in the dark in a rolling hybridization incubator.
NOTE: Amplification time can be increased beyond 4 h if desired.
NOTE: Amplification may alternatively be carried out at room temperature.
6. Remove excess hairpins by washing in pre-heated $5 \times$ SSCT (volume $2V$ per wash): 2×15 min at 37 °C in the dark.
NOTE: If amplification is performed at room temperature, hairpin washes should also be performed at room temperature.

Buffer recipes for sample preparation

10% SDS

10% SDS

For 100 mL of solution

10 g of SDS powder

Fill up to 100 mL with ultrapure H₂O

Heat solution to 68 °C to dissolve powder if necessary.

2× SSC, 0.1% SDS

2× sodium chloride sodium citrate (SSC)

0.1% SDS

For 100 mL of solution

10 mL of 20× SSC

1 mL of 10% SDS

Fill up to 100 mL with ultrapure H₂O

0.2× SSC, 0.1% SDS

0.2× sodium chloride sodium citrate (SSC)

0.1% SDS

For 100 mL of solution

1 mL of 20× SSC

1 mL of 10% SDS

Fill up to 100 mL with ultrapure H₂O

5× SSCT

5× sodium chloride sodium citrate (SSC)

0.1% Tween 20

For 40 mL of solution

10 mL of 20× SSC

400 µL of 10% Tween 20

Fill up to 40 mL with ultrapure H₂O

Reagents and supplies

NorthernMax™ denaturing gel buffer (10X) (ThermoFisher Scientific, Cat. # AM8676)
NorthernMax™ (10X) running buffer (ThermoFisher Scientific, Cat. # AM8671)
Nylon membrane, positively-charged (Roche, Cat. # 11209272001)
3 mm Whatman chromatography paper (GE Healthcare Life Sciences, Cat. # 3030-6188)
Whatman TurboBlotter (Fisher Scientific, Cat. # 09-301-184)
Whatman gel blotting paper (GE Healthcare Life Sciences, Cat. # 10427805)
Whatman GB003 filter paper (GE Healthcare Life Sciences, Cat. # 10427812)
Hybridization bottle (Wheaton, Cat. # 805000 or 805021)
Rolling hybridization incubator (SciGene, Cat. # 1040-50-1)

S1 HCR™ Technology Citation Notes

For citation, please select from the list below as appropriate for your application:

- **HCR™ RNA-ISH**

HCR™ RNA in situ hybridization (RNA-ISH) offers unmatched performance, robustness, and versatility imaging RNA targets in diverse organisms and sample types (Choi et al., 2010, Choi et al., 2014, Choi et al., 2018):

- **HCR™ RNA-FISH**

HCR™ RNA-FISH enables 10-plex, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) with automatic background suppression throughout the protocol for dramatically enhanced performance (signal-to-background, quantitative precision, single-molecule fidelity) and ease-of-use (no probe set optimization for new targets and organisms).

- **Enzymatic HCR™ RNA-CISH/RNA-FISH**

Enzymatic HCR™ RNA-ISH integrates enzymatic signal amplification to enable extreme-sensitivity RNA imaging using either chromogenic or fluorescent staining (RNA-CISH/RNA-FISH). In tissue sections, entirely protease-free workflows preserve sample morphology and maintain protein target integrity, enabling seamless compatibility with existing immunohistochemistry (IHC)/immunofluorescence (IF) assays. HCR™ RNA-CISH offers the convenience of brightfield microscopy and the option of archival staining.

- **10-Plex HCR™ Spectral Imaging**

HCR™ RNA-FISH/IF enables quantitative high-resolution imaging of 10 RNA and/or protein targets with 1-step HCR™ signal amplification for all targets simultaneously. The method is suitable even for whole-mounts and delicate samples as it requires no repeated staining, imaging, registration, or stripping (Schulte et al., 2024).

- **HCR™ RNA-FISH/IF**

HCR™ RNA-FISH/IF enables a unified approach to multiplex, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) and protein immunofluorescence (IF), with quantitative 1-step enzyme-free signal amplification performed for all RNA and protein targets simultaneously (Schwarzkopf et al., 2021).

- **HCR™ IF**

HCR™ IF enables multiplex, quantitative, high-resolution protein immunofluorescence (IF) in highly autofluorescent samples (e.g., FFPE brain tissue sections) (Schwarzkopf et al., 2021).

- **Subcellular Quantitative RNA and Protein Imaging**

HCR™ RNA-FISH enables analog relative quantitation of RNA and/or protein targets with subcellular resolution in the anatomical context of thick autofluorescent samples (e.g., whole-mount vertebrate embryos) (Trivedi et al., 2018, Choi et al., 2018, Schwarzkopf et al., 2021).

- **Single-Molecule Quantitative RNA Imaging**

HCR™ RNA-FISH enables digital RNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) (Shah et al., 2016, Choi et al., 2018).

- **Read-Out/Read-In Analysis Framework**

The read-out/read-in analysis framework enables bidirectional quantitative discovery in an anatomical context (Trivedi et al., 2018).

- **Protocols in Diverse Sample Types**

Protocols for HCR™ RNA-FISH and/or IF in diverse sample types are adapted from the zoo paper ([Choi et al., 2016](#)):

- bacteria in suspension
- FFPE human tissue sections
- generic sample in solution
- generic sample on a slide
- mammalian cells on a slide
- mammalian cells in suspension
- whole-mount chicken embryos
- whole-mount fruit fly embryos
- whole-mount mouse embryos
- whole-mount nematode larvae
- whole-mount sea urchin embryos
- whole-mount zebrafish embryos and larvae

- **HCR™ RNA Flow Cytometry**

HCR™ RNA Flow Cytometry enables analog RNA relative quantitation for high-throughput expression profiling of mammalian cells and bacteria without the need to engineer reporter lines ([Choi et al., 2018](#)).

- **HCR™ Northern Blots**

HCR™ Northern Blots enable simultaneous quantification of RNA target size and abundance with automatic background suppression throughout the protocol ([Schwarzkopf & Pierce, 2016](#)).

- **HCR™ Amplifiers**

HCR™ Amplifiers enable multiplex, quantitative, 1-step, isothermal, enzyme-free signal amplification in diverse technological settings ([Dirks & Pierce, 2004](#)).