

## HCR™ RNA-FISH protocol for whole-mount nematode larvae (*Caenorhabditis elegans*)

This protocol has not been optimized for all stages and should only be used as a template.

### Technical support



[support@molecularinstruments.com](mailto:support@molecularinstruments.com)

### Safety data sheets (SDS)

[www.molecularinstruments.com/safety](http://www.molecularinstruments.com/safety)

### Patents

Molecular Instruments® products are protected by and for use under patents: [www.molecularinstruments.com/patents](http://www.molecularinstruments.com/patents)



## **HCR™ RNA-FISH**

*Multiplexed, quantitative, high-resolution RNA imaging*

### **Multiplexed Experiment**

- Order one HCR™ RNA-FISH bundle per target RNA

### **Example 2-Plex Experiment**

- HCR™ RNA-FISH bundle for target mRNA1
  - HCR™ probe set: target mRNA1 for use with amplifier B1
  - HCR™ amplifier: B1-647
  - HCR™ RNA-FISH buffers: probe hybridization buffer, probe wash buffer, amplification buffer (for use with all bundles)
- HCR™ RNA-FISH bundle for target mRNA2
  - HCR™ probe set: target mRNA2 for use with amplifier B2
  - HCR™ amplifier: B2-488

### **Storage conditions**

- Store HCR™ probe sets, HCR™ amplifiers, HCR™ RNA-FISH probe hybridization buffer, and HCR™ RNA-FISH probe wash buffer at -20 °C.
- Store HCR™ amplification buffer at 4 °C.
- On the bench top, keep stock solutions on ice.
- Make sure all solutions are well mixed before use.

## Preparation of fixed whole-mount nematode larvae

1. Wash nematode larvae off plate with 1 mL of M9 buffer and transfer to a 1.5 mL tube.
2. Centrifuge at  $200 \times g$  for 2 min to bring larvae to the bottom of the tube and remove the M9 buffer.
3. Wash larvae 3 times with 1 mL of M9 buffer each. Centrifuge at  $200 \times g$  for 2 min between washes.
4. Centrifuge and remove  $\approx 800 \mu\text{L}$  of M9 buffer.
5. Aliquot larvae sufficient for a reaction into 1.5 mL tubes.
6. Add 1 mL of 4% paraformaldehyde (PFA).  
*CAUTION: use PFA with extreme care as it is a hazardous material.*  
*NOTE: use fresh PFA to avoid increased autofluorescence.*
7. Immediately freeze sample at  $-80^\circ\text{C}$  overnight before use.  
*NOTE: Larvae could stay in  $-80^\circ\text{C}$  freezer for long-term storage.*
8. Fix larvae by thawing at room temperature for 45 min.
9. Wash larvae 2 times with 1 mL of PBST each. Bring larvae to the bottom of the tube with centrifugation at  $200 \times g$  for 2 min in between washes.
10. Treat larvae with 1 mL of proteinase K ( $100 \mu\text{g}/\text{mL}$ ) for 10 min at  $37^\circ\text{C}$ .  
*NOTE: Proteinase K concentration and treatment time should be reoptimized for each batch of proteinase K, or for samples at a different developmental stage.*
11. Wash larvae 2 times with 1 mL of PBST each.
12. Incubate in 1 mL of 2 mg/mL of glycine solution for 15 min on ice.
13. Wash larvae 2 times with 1 mL of PBST each.

## Multiplexed HCR™ RNA-FISH protocol

### Detection stage

1. Incubate larvae in 1 mL of 50% PBST / 50% probe hybridization buffer for 5 min at room temperature.  
*CAUTION: probe hybridization buffer contains formamide, a hazardous material.*
2. Centrifuge at  $200 \times g$  for 2 min to remove solution.
3. Pre-hybridize larvae in 300  $\mu\text{L}$  of probe hybridization buffer at 37 °C for 1 h.
4. Prepare probe solution by adding 2 pmol of each probe set (e.g. 2  $\mu\text{L}$  of 1  $\mu\text{M}$  stock) to 200  $\mu\text{L}$  of probe hybridization buffer at 37 °C.  
*NOTE: For single-molecule RNA imaging, use higher probe concentration (e.g., 16 nM) to increase probe hybridization yield. If desired, this approach can also be used to increase signal for subcellular quantitative RNA imaging.*
5. Add the probe solution to reach a final hybridization volume of 500  $\mu\text{L}$ .
6. Incubate larvae overnight (>12 h) at 37 °C.
7. Remove excess probes by washing larvae  $4 \times 15$  min with 1 mL of probe wash buffer at 37 °C.  
*CAUTION: probe wash buffer contains formamide, a hazardous material.*  
*NOTE: pre-heat probe wash buffer to 37 °C before use.*  
*NOTE: bring larvae to the bottom of the tube with centrifugation at  $500 \times g$  for 2 min for each wash.*
8. Wash larvae  $2 \times 5$  min with 1 mL of  $5 \times$  SSCT at room temperature.

### Amplification stage

1. Pre-amplify larvae with 300  $\mu\text{L}$  of amplification buffer for 30 min at room temperature.  
*NOTE: equilibrate amplification buffer to room temperature before use.*
2. Separately prepare 30 pmol of hairpin h1 and 30 pmol of hairpin h2 by snap cooling 10  $\mu\text{L}$  of 3  $\mu\text{M}$  stock (heat at 95 °C for 90 seconds and cool to room temperature in a dark drawer for 30 min).  
*NOTE: Hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. h1 and h2 should be snap cooled in separate tubes.*
3. Prepare hairpin solution by adding snap-cooled h1 hairpins and snap-cooled h2 hairpins to 200  $\mu\text{L}$  of amplification buffer at room temperature.
4. Add the hairpin solution to reach a final amplification volume of 500  $\mu\text{L}$ .
5. Incubate the larvae overnight (>12 h) in the dark at room temperature.  
*NOTE: For single-molecule RNA imaging, amplify for a shorter period of time to ensure single-molecule dots are diffraction-limited.*
6. Remove excess hairpins by washing with 1 mL of  $5 \times$  SSCT at room temperature:
  - (a)  $2 \times 5$  min
  - (b)  $2 \times 30$  min

(c)  $1 \times 5$  min

*NOTE: bring larvae to the bottom of the tube with centrifugation at  $500 \times g$  for 2 min for each wash.*

7. Samples can be stored at  $4^\circ\text{C}$  protected from light before microscopy.

## Buffer recipes

### M9 buffer

22 mM  $\text{KH}_2\text{PO}_4$   
42 mM  $\text{Na}_2\text{HPO}_4$   
20.5 mM NaCl  
1 mM  $\text{MgSO}_4$

### For 1 L of solution

3 g of  $\text{KH}_2\text{PO}_4$   
6 g of  $\text{Na}_2\text{HPO}_4$   
5 g of NaCl  
1 mL of 1 M  $\text{MgSO}_4$   
Fill up to 1 L with ultrapure  $\text{H}_2\text{O}$   
Sterilize by autoclaving  
Store buffer at 4 °C before use

### 4% Paraformaldehyde (PFA)

4% PFA  
1× PBS

### For 40 mL of solution

10 mL of 16% PFA solution  
4 mL of 10× PBS  
Fill up to 40 mL with ultrapure  $\text{H}_2\text{O}$

### PBST

1× PBS  
0.1% Tween 20

### For 50 mL of solution

5 mL of 10× PBS  
500  $\mu\text{L}$  of 10% Tween 20  
Fill up to 50 mL with ultrapure  $\text{H}_2\text{O}$

### Proteinase K solution

100  $\mu\text{g}/\text{mL}$  proteinase K

### For 1 mL of solution

5  $\mu\text{L}$  of 20 mg/mL proteinase K  
Fill up to 1 mL with PBST

### Glycine solution

2 mg/mL glycine  
PBST

### For 50 mL of solution

100 mg of glycine  
Fill up to 50 mL with PBST

### 5× SSCT

5× sodium chloride sodium citrate (SSC)  
0.1% Tween 20

### For 40 mL of solution

10 mL of 20× SSC  
400  $\mu\text{L}$  of 10% Tween 20  
Fill up to 40 mL with ultrapure  $\text{H}_2\text{O}$

*NOTE: avoid using calcium chloride and magnesium chloride in PBS as this leads to increased autofluorescence in the samples.*

## HCR™ Technology Citation Notes

For citation, please select from the list below as appropriate for your application:

- **HCR™ IF + HCR™ RNA-FISH**

[HCR™ IF + HCR™ RNA-FISH](#) enables a unified approach to multiplexed, quantitative, high-resolution protein immunofluorescence (IF) and RNA fluorescence in situ hybridization (RNA-FISH), with quantitative 1-step enzyme-free signal amplification performed for all protein and RNA targets simultaneously ([Schwarzkopf et al., 2021](#)).

- **HCR™ IF**

[HCR™ IF](#) enables multiplexed, quantitative, high-resolution protein immunofluorescence (IF) in highly autofluorescent samples (e.g., FFPE brain tissue sections) ([Schwarzkopf et al., 2021](#)).

- **HCR™ RNA-FISH (v3.0)**

Third-generation [HCR™ RNA-FISH \(v3.0\)](#) enables multiplexed, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) with [automatic background suppression throughout the protocol](#) for dramatically enhanced performance (signal-to-background, subcellular quantitative RNA imaging precision, single-molecule quantitative RNA imaging fidelity) and ease-of-use (no probe set optimization for new targets and organisms) ([Choi et al., 2018](#)). Quantitative analysis modes:

- [Subcellular quantitative RNA imaging](#): analog mRNA relative quantitation with subcellular resolution in the anatomical context of thick autofluorescent samples.
- [Single-molecule quantitative RNA imaging](#): digital mRNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples.
- [Quantitative RNA flow cytometry](#): analog mRNA relative quantitation for high-throughput expression profiling of mammalian cells and bacteria.

[Protocols for HCR™ RNA-FISH \(v3.0\)](#) in diverse organisms are adapted from the Zoo paper.

- **Subcellular quantitative RNA imaging**

[Subcellular quantitative RNA imaging](#) enables mRNA relative quantitation with subcellular resolution in the anatomical context of thick autofluorescent samples (e.g., whole-mount vertebrate embryos). The [read-out/read-in analysis framework](#) enables bidirectional quantitative discovery in an anatomical context ([Trivedi et al., 2018](#)).

- **Zoo paper**

Protocols for multiplexed mRNA imaging in diverse sample types ([Choi et al., 2016](#)):

- bacteria in suspension
- FFPE human tissue sections
- generic sample in solution
- generic sample on a slide
- mammalian cells on a slide
- mammalian cells in suspension
- whole-mount chicken embryos
- whole-mount fruit fly embryos
- whole-mount mouse embryos
- whole-mount nematode larvae
- whole-mount sea urchin embryos
- whole-mount zebrafish embryos and larvae

- **Single-molecule quantitative RNA imaging**

[Single-molecule quantitative RNA imaging](#) enables RNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) ([Shah et al., 2016](#)).

- **HCR™ northern blots**

[HCR™ northern blots](#) enable multiplexed quantification of RNA target size and abundance for up to 5 target RNAs ([Schwarzkopf & Pierce, 2016](#)).

- **HCR™ RNA-FISH (v2.0)**

Second-generation in situ HCR™ RNA-FISH technology (v2.0) using DNA HCR™ probes and DNA HCR™ amplifiers: 10× increase in signal, 10× reduction in cost, dramatic increase in reagent durability ([Choi et al., 2014](#)).

- **HCR™ RNA-FISH (v1.0)**

First-generation HCR™ RNA-FISH technology (v1.0) using RNA HCR™ probes and RNA HCR™ amplifiers: multiplexed mRNA imaging in whole-mount vertebrate embryos with simultaneous signal amplification for up to 5 target mRNAs ([Choi et al., 2010](#)).

- **HCR™ technology**

[HCR™ amplifiers](#) enable multiplexed, quantitative, 1-step, isothermal, enzyme-free signal amplification in diverse technological settings ([Dirks & Pierce, 2004](#)).