

## HCR™ RNA flow cytometry protocol for bacteria in suspension

This protocol has been optimized for *Escherichia coli* and should only be used as a template for other types of bacteria.

### Technical support



[support@molecularinstruments.com](mailto:support@molecularinstruments.com)

### Safety data sheets (SDS)

[www.molecularinstruments.com/safety](http://www.molecularinstruments.com/safety)

### Patents

Molecular Instruments® products are protected by and for use under patents: [www.molecularinstruments.com/patents](http://www.molecularinstruments.com/patents)



## HCR™ RNA Flow Cytometry

*Multiplexed quantitative RNA flow cytometry*

### Multiplexed Experiment

- Order one HCR™ RNA-FISH bundle per target RNA

### Example 2-Plex Experiment

- HCR™ RNA-FISH bundle for target mRNA1
  - HCR™ probe set: target mRNA1 for use with amplifier B1
  - HCR™ amplifier: B1-647
  - HCR™ RNA-FISH buffers: probe hybridization buffer, probe wash buffer, amplification buffer (for use with all bundles)
- HCR™ RNA-FISH bundle for target mRNA2
  - HCR™ probe set: target mRNA2 for use with amplifier B2
  - HCR™ amplifier: B2-488

### Storage conditions

- Store HCR™ probe sets, HCR™ amplifiers, HCR™ RNA-FISH probe hybridization buffer, and HCR™ RNA-FISH probe wash buffer at -20 °C.
- Store HCR™ amplification buffer at 4 °C.
- On the bench top, keep stock solutions on ice.
- Make sure all solutions are well mixed before use.

## Sample preparation

1. Grow *E. coli* from streaked plate or frozen glycerol stocks in 2–3 mL of LB media overnight in a 37 °C shaker.
2. Dilute to make a 5 mL liquid culture with  $OD_{600} = 0.05$ .
3. Incubate in a 37 °C shaker until  $OD_{600} \approx 0.5$  (exponential phase).
4. Aliquot 1 mL of cells and centrifuge for 10 min.  
*NOTE: Centrifugation should be as gentle as possible to pellet cells. For E. coli, all centrifugation steps are done at  $4000 \times g$ .*
5. Remove supernatant and re-suspend cells in 750  $\mu\text{L}$  of  $1 \times$  phosphate-buffered saline (PBS).  
*NOTE: remove all solutions via pipetting throughout the protocol.*
6. Add 250  $\mu\text{L}$  of 4% formaldehyde and incubate overnight at 4 °C.  
*CAUTION: use formaldehyde with extreme care as it is a hazardous material.*
7. Centrifuge for 10 min and remove supernatant.
8. Re-suspend cells in 150  $\mu\text{L}$  of  $1 \times$  PBS.
9. Add 850  $\mu\text{L}$  of 100% MeOH and store cells at -20 °C before use.

## Multiplexed HCR™ RNA flow cytometry protocol

### Detection stage

1. Transfer 150  $\mu\text{L}$  of cells into a 1.5 mL eppendorf tube.
2. Centrifuge for 5 min and remove supernatant.
3. Wash cells with 500  $\mu\text{L}$  of 1 $\times$  PBST. Centrifuge for 5 min and remove the supernatant.
4. Re-suspend the pellet with 400  $\mu\text{L}$  of probe hybridization buffer and pre-hybridize for 1 h at 37 °C.  
*CAUTION: probe hybridization buffer contains formamide, a hazardous material.*  
*NOTE: Probe hybridization buffer for cells in suspension will be used in this protocol.*
5. In the meantime, prepare probe solution by adding 2 pmol of each probe set (e.g. 2  $\mu\text{L}$  of 1  $\mu\text{M}$  stock) to 100  $\mu\text{L}$  of probe hybridization buffer pre-heated to 37 °C.
6. Add the probe solution directly to the sample to reach a final probe concentration of 4 nM.  
*NOTE: use higher probe concentration (e.g., 16 nM) to improve signal for quantitative RNA flow cytometry.*
7. Incubate the sample overnight (>12 h) at 37 °C.
8. Add 1 mL of probe wash buffer to the sample to dilute the probes.  
*CAUTION: probe wash buffer contains formamide, a hazardous material.*  
*NOTE: pre-heat probe wash buffer to 37 °C before use.*
9. Centrifuge for 5 min and remove the wash solution.
10. Re-suspend the cell pellet with 500  $\mu\text{L}$  of probe wash buffer (pre-heated to 37 °C).
11. Incubate for 5 min at 37 °C and remove the wash solution by centrifugation for 5 min.
12. Repeat steps 10 and 11 for two additional times but with 10 min incubation.
13. Proceed to hairpin amplification.

### Amplification stage

1. Re-suspend the cell pellet with 150  $\mu\text{L}$  of amplification buffer and pre-amplify for 30 min at room temperature.  
*NOTE: Amplification buffer for cells in suspension will be used in this protocol. Equilibrate amplification buffer to room temperature before use.*
2. Separately prepare 15 pmol of hairpin h1 and 15 pmol of hairpin h2 by snap cooling 5  $\mu\text{L}$  of 3  $\mu\text{M}$  stock (heat at 95 °C for 90 seconds and cool to room temperature in a dark drawer for 30 min).  
*NOTE: Hairpins h1 and h2 are provided in hairpin storage buffer ready for snap cooling. h1 and h2 should be snap cooled in separate tubes.*
3. Prepare hairpin solution by adding snap-cooled h1 hairpins and snap-cooled h2 hairpins to 100  $\mu\text{L}$  of amplification buffer at room temperature.
4. Add the hairpin solution directly to the sample to reach a final hairpin concentration of 60 nM.
5. Incubate the sample overnight (>12 h) in the dark at room temperature.

6. Add 1 mL of 5× SSCT at room temperature to the sample to dilute the hairpins.
7. Centrifuge for 5 min and remove the wash solution.
8. Re-suspend the cell pellet with 500  $\mu$ L of 5× SSCT.
9. Incubate for 5 min at room temperature and remove the wash solution by centrifugation for 5 min.
10. Repeat steps 8 and 9 for two additional times but with 10 min incubation.
11. Re-suspend the cell pellet in desired buffer (e.g. 5× SSCT) and volume.
12. Samples can be stored at 4 °C protected from light before microscopy or flow cytometry.

*NOTE: filter cells before flow cytometry.*

## Buffer recipes

### **LB media**

5 g of Novagen LB Broth Miller powder  
Fill up to 200 mL with ultrapure H<sub>2</sub>O  
Autoclave at 121 °C for 20 min

### **4% formaldehyde**

4% formaldehyde  
1 × PBS

### **For 4 mL of solution**

1 mL of 16% formaldehyde  
0.4 mL of 10 × PBS  
Fill up to 4 mL with ultrapure H<sub>2</sub>O

### **5 × SSCT**

5 × sodium chloride sodium citrate (SSC)  
0.1% Tween 20

### **For 40 mL of solution**

10 mL of 20 × SSC  
400 μL of 10% Tween 20  
fill up to 40 mL with ultrapure H<sub>2</sub>O

*NOTE: avoid using calcium chloride and magnesium chloride in PBS as this leads to increased autofluorescence in the samples.*

## HCR™ Technology Citation Notes

For citation, please select from the list below as appropriate for your application:

- **HCR™ IF + HCR™ RNA-FISH**

[HCR™ IF + HCR™ RNA-FISH](#) enables a unified approach to multiplexed, quantitative, high-resolution protein immunofluorescence (IF) and RNA fluorescence in situ hybridization (RNA-FISH), with quantitative 1-step enzyme-free signal amplification performed for all protein and RNA targets simultaneously ([Schwarzkopf et al., 2021](#)).

- **HCR™ IF**

[HCR™ IF](#) enables multiplexed, quantitative, high-resolution protein immunofluorescence (IF) in highly autofluorescent samples (e.g., FFPE brain tissue sections) ([Schwarzkopf et al., 2021](#)).

- **HCR™ RNA-FISH (v3.0)**

Third-generation [HCR™ RNA-FISH \(v3.0\)](#) enables multiplexed, quantitative, high-resolution RNA fluorescence in situ hybridization (RNA-FISH) with [automatic background suppression throughout the protocol](#) for dramatically enhanced performance (signal-to-background, subcellular quantitative RNA imaging precision, single-molecule quantitative RNA imaging fidelity) and ease-of-use (no probe set optimization for new targets and organisms) ([Choi et al., 2018](#)). Quantitative analysis modes:

- [Subcellular quantitative RNA imaging](#): analog mRNA relative quantitation with subcellular resolution in the anatomical context of thick autofluorescent samples.
- [Single-molecule quantitative RNA imaging](#): digital mRNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples.
- [Quantitative RNA flow cytometry](#): analog mRNA relative quantitation for high-throughput expression profiling of mammalian cells and bacteria.

[Protocols for HCR™ RNA-FISH \(v3.0\)](#) in diverse organisms are adapted from the Zoo paper.

- **Subcellular quantitative RNA imaging**

[Subcellular quantitative RNA imaging](#) enables mRNA relative quantitation with subcellular resolution in the anatomical context of thick autofluorescent samples (e.g., whole-mount vertebrate embryos). The [read-out/read-in analysis framework](#) enables bidirectional quantitative discovery in an anatomical context ([Trivedi et al., 2018](#)).

- **Zoo paper**

Protocols for multiplexed mRNA imaging in diverse sample types ([Choi et al., 2016](#)):

- bacteria in suspension
- FFPE human tissue sections
- generic sample in solution
- generic sample on a slide
- mammalian cells on a slide
- mammalian cells in suspension
- whole-mount chicken embryos
- whole-mount fruit fly embryos
- whole-mount mouse embryos
- whole-mount nematode larvae
- whole-mount sea urchin embryos
- whole-mount zebrafish embryos and larvae

- **Single-molecule quantitative RNA imaging**

[Single-molecule quantitative RNA imaging](#) enables RNA absolute quantitation with single-molecule resolution in the anatomical context of thick autofluorescent samples (e.g., 0.5 mm adult mouse brain sections) ([Shah et al., 2016](#)).

- **HCR™ northern blots**

[HCR™ northern blots](#) enable multiplexed quantification of RNA target size and abundance for up to 5 target RNAs ([Schwarzkopf & Pierce, 2016](#)).

- **HCR™ RNA-FISH (v2.0)**

Second-generation in situ HCR™ RNA-FISH technology (v2.0) using DNA HCR™ probes and DNA HCR™ amplifiers: 10× increase in signal, 10× reduction in cost, dramatic increase in reagent durability ([Choi et al., 2014](#)).

- **HCR™ RNA-FISH (v1.0)**

First-generation HCR™ RNA-FISH technology (v1.0) using RNA HCR™ probes and RNA HCR™ amplifiers: multiplexed mRNA imaging in whole-mount vertebrate embryos with simultaneous signal amplification for up to 5 target mRNAs ([Choi et al., 2010](#)).

- **HCR™ technology**

[HCR™ amplifiers](#) enable multiplexed, quantitative, 1-step, isothermal, enzyme-free signal amplification in diverse technological settings ([Dirks & Pierce, 2004](#)).